



**Single-Dose Packaging  
Solution with Real-Time  
Monitoring to  
Boost Patient  
Adherence**



# Boosting patient adherence without adding to patient burden

Innovative CubixxDT transforms adherence issues into cost-effective exception management, delivering value in a wholly disposable structure.


## How CubixxDT Works



1 Provider receives request and prepares dose pack



2 Apply circuit label, attach transmitter and dispense to patient




Patient 3B

3 Patient removes dose



4 Cellular information transmitted in real-time to defined stakeholders



5 Customized portal reports on patients and participants

# Delivers Value Without Disrupting Patients or Providers




Self-contained, cost-effective and a wholly disposable infrastructure in each pack



No gadget return



No network setup



No additional structure

---

## Improves Compliance Without Requiring Additional Actions from Patients

- » Transforms adherence issues into cost-efficient exception management
- » Disposable, cellular connectivity means no patient record keeping, diary, app download or cell-phone required
- » Patient-centric innovation built around patient behaviors
- » Records and communicates adherence information in real-time

---

## Provides Immediate Visibility into HIPAA-Compliant Data to Improve Care Planning

- » No wait for utilization data
- » Online, centralized dashboard for oversight of one or multiple locations
- » Cellular information transmitted in real-time to defined stakeholders
- » Customized portal and reports generation on patients and participants

## About Us

Applying everything from RFID innovation to remote monitoring tools to customized user interfaces that capture real-time patient data, we're transforming the way healthcare providers, pharmaceutical manufacturers, payers and patients approach medication storage and dispensing.

For additional information on Cubixx Solutions, please visit [cubixxsolutions.com](http://cubixxsolutions.com) or call **844.428.2499**, extension 7906.